

For Research Use Only



# MultiPro™ 5CFLX Mouse IgG1 Isotype Control (1F8D3)

Catalog Number: G66360-1-5C

## Basic Information

Catalog Number:	GenBank Accession Number:	CloneNo.:
G66360-1-5C	GeneID (NCBI):	1F8D3
Size:	Full Name:	Conjugate:
500ug/mL	MultiPro™ 5CFLX Mouse IgG1 Isotype Control (1F8D3)	5CFLX
Source:		Full Oligo Sequence:
Mouse		CGGAGATGTGTATAAGACAGTCTG AATATGTCGTGCCCATATAAGAAA
Isotype:		Barcode Sequence:
IgG1		TCTGAATATGTCGTG

## Applications

Tested Applications:  
Single Cell (Intra)  
Species Specificity:  
NA

## Background Information

This antibody is used as the isotype control of mouse IgG1 (Kappa). It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

## Storage

Storage:  
2-8°C  
Storage Buffer:  
PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



G66360-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G66360-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.