

MultiPro[™] 5CFLX Mouse IgG2a Isotype Control (11A1B2)

Catalog Number: **G66360-2-5C**

Basic Information

Catalog Number:

G66360-2-5C

Size:

500ug/mL

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

Full Name:

MultiPro[™] 5CFLX Mouse IgG2a

Isotype Control (11A1B2)

CloneNo.:

11A1B2

Conjugate:

5CFLX

Full Oligo Sequence:

CGGAGATGTGTATAAGAGACAGAGAT

ACTCGTGTACCCCATATAAGAAA

Barcode Sequence:

AGATACTCGTGTAC

Applications

Tested Applications:

Single Cell (Intra)

Species Specificity:

NA

Background Information

This antibody is used as the isotype control of mouse IgG2a. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Storage


Storage:

2-8°C


Storage Buffer:

PBS with 1mM EDTA and 0.09% sodium azide

Selected Validation Data



G66360-2-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G66360-2-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.