

For Research Use Only



MultiPro™ 5CFLX Mouse IgG3 Isotype Control (1H4A5)

Catalog Number: G66360-4-5C

Basic Information

Catalog Number:

G66360-4-5C

Size:

500ug/mL

Source:

Mouse

Isotype:

IgG3

GenBank Accession Number:

GeneID (NCBI):

Full Name:

MultiPro™ 5CFLX Mouse IgG3 Isotype Control (1H4A5)

CloneNo.:

1H4A5

Conjugate:

5CFLX

Full Oligo Sequence:

CGGAGATGTGTATAAGAGACAGCTCC
TACAGGACAGACCCATATAAGAAA

Barcode Sequence:

CTCCTACAGGACAGA

Applications

Tested Applications:

Single Cell (Intra)

Species Specificity:

NA

Background Information

This antibody is used as the isotype control of mouse IgG3. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Storage

Storage:

2-8°C

Storage Buffer:

PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G66360-4-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G66360-4-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.