

For Research Use Only

MultiPro™ 5CFLX Anti-Human TGM2 (2D4C11)



Catalog Number: G68006-1-5C

Basic Information

Catalog Number:	GenBank Accession Number:	CloneNo.:
G68006-1-5C	BC003551	2D4C11
Size:	GeneID (NCBI):	Conjugate:
500ug/mL	7052	5CFLX
Source:	ENSEMBL Gene ID:	Full Oligo Sequence:
Mouse	ENSG00000198959	CGGAGATGTGTATAAGACAGTATA
Isotype:	UNIPROT ID:	CTGAAGGATATCCCATATAAGAAA
IgG2a	P21980	Barcode Sequence:
Immunogen Catalog Number:	Full Name:	TATACTGAAGGATAT
AG7462	MultiPro™ 5CFLX Anti-Human TGM2 (2D4C11)	

Applications

Tested Applications:
Single Cell (Intra)
Species Specificity:
Human

Background Information

Storage

Storage:
2-8°C
Storage Buffer:
PBS with 1mM EDTA and 0.09% sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



G68006-1-5C was used to stain Resting PBMC and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.



G68006-1-5C was used to stain PBMC under 4hr LPS + TI treatment and analyzed with 10x Genomics Gene Expression Flex with Feature Barcodes and Multiplexing kit with Fix-Stain protocol.