

For Research Use Only

HRP-conjugated IgG Lambda Light Chain Polyclonal antibody



Catalog Number: **HRP-11541**

Basic Information

Catalog Number:

HRP-11541

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2064

GenBank Accession Number:

BC020233

GeneID (NCBI):

3535

Full Name:

immunoglobulin lambda locus

Calculated MW:

237 aa, 25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB

Species Specificity:

human

Positive Controls:

WB : human plasma,

Background Information

This antibody detects the lambda light chain of human IgG.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human plasma were subjected to SDS PAGE followed by western blot with HRP-11541 (IgG Lambda Light Chain antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.