

For Research Use Only

# HRP-conjugated IgA Monoclonal antibody



Catalog Number: HRP-60099

## Basic Information

<b>Catalog Number:</b> HRP-60099	<b>GenBank Accession Number:</b> BC016369	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 3493	<b>CloneNo.:</b> 3A5A5
<b>Source:</b> Mouse	<b>Full Name:</b> immunoglobulin heavy constant alpha 1	<b>Recommended Dilutions:</b> WB 1:5000-1:50000
<b>Isotype:</b> IgG1	<b>Calculated MW:</b> 496 aa, 53 kDa	
<b>Immunogen Catalog Number:</b> AG7417	<b>Observed MW:</b> 55-60 kDa	

## Applications

<b>Tested Applications:</b> WB	<b>Positive Controls:</b> WB : human saliva, human placenta tissue, human plasma
<b>Species Specificity:</b> human	

## Background Information

Ig alpha is the major immunoglobulin class in body secretions. It may serve both to defend against local infection and to prevent access of foreign antigens to the general immunologic system. This antibody recognizes the Chain C region of human Ig alpha-1 (IGHA1).

## Storage

**Storage:**  
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
**Aliquoting is unnecessary for -20°C storage**

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



50 ug, 25 ug, 12.5 ug human saliva were subjected to SDS PAGE followed by western blot with HRP-60099 (Human IgA antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.