

For Research Use Only

HRP-conjugated Human IgG Kappa chain Monoclonal antibody



Catalog Number:HRP-67761

1 Publications

Basic Information

Catalog Number:

HRP-67761

Size:

1000 μ g/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

155 kDa

Observed MW:

~25 kDa

Purification Method:

Protein G purification

CloneNo.:

4D11F7

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB

Cited Applications:

WB

Species Specificity:

Human

Positive Controls:

WB : human plasma,

Background Information

This antibody detects the Kappa chain of human IgG,

Notable Publications

Author	Pubmed ID	Journal	Application
Ying Liu	38347502	Pharm Biol	WB

Storage

Storage:

Store at -20°C. Avoid exposure to light.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with HRP-67761 (Human IgG Kappa chain antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.