

For Research Use Only

HRP-conjugated IgG lambda chain Monoclonal antibody



Catalog Number: HRP-67762

Basic Information

Catalog Number: HRP-67762	GenBank Accession Number: BC007782	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 3537	CloneNo.: 1D12B5
Source: Mouse	Full Name: immunoglobulin lambda constant 1 (Mcg marker)	Recommended Dilutions: WB 1:2000-1:10000
Isotype: IgG2a	Calculated MW: 233 aa, 25 kDa	
Immunogen Catalog Number: AG14783	Observed MW: ~25 kDa	

Applications

Tested Applications: WB	Positive Controls: WB : human plasma,
Species Specificity: Human	

Background Information

This antibody detects Lambda chain of human IgG.

Storage

Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with HRP-67762 (Human IgG lambda chain antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.