

Catalog Number: HZ-1323

Data Sheet

HumanKine® IL-15 (Recombinant Human)



Animal Component-Free

Human cell expressed

Tag-Free

Endotoxin Free


Product Description

Alternative Names	Interleukin 15, IL15, IL-15MGC9721, Interleukin-15
Source	Human Embryonic Kidney cells (HEK293). HEK293-derived IL-15 protein


Specifications

Test	Method	Specification
Activity	Dose-dependent proliferation of the NK-92 human natural killer cell line.	0.07-0.37 ng/mL for NK-92 bioassay.
Molecular Mass	SDS-PAGE	12 to 24 kDa reduced, 11 to 23 kDa non-reduced, monomer, glycosylated
Purity	SDS-PAGE	>95%
Endotoxin	LAL	<1 EU/μg

Activity Data



SDS-PAGE



www.ptglab.com

Document #: FR-QA118-101 Rev 0
Data Sheet Version #: 1

Proteintech Group, Inc.

5500 Pearl Street, Suite 400 Rosemont, IL 60012
t: 1-888-478-4522; f: 1-312-455-8408
Email: proteintech@ptglab.com

Preparation	
Shipping Temperature	ambient temperature
Formulation	50 mM Acetate pH 4.0 + 150mM NaCl
Reconstitution	Briefly centrifuge the vial before opening. It is recommended to reconstitute the protein to 0.2 mg/mL in sterile 1x PBS pH 7.4 containing 0.1% endotoxin-free recombinant human serum albumin (HSA). Gently swirl or tap vial to mix.

Stability and Storage	Product Form	Temperature Conditions	Storage Time (From Date of Receipt)
	Lyophilized	-20°C to -80°C	Until Expiry Date
	Lyophilized	Room Temperature	2 weeks
	Reconstituted as per CofA	-20°C to -80°C	6 months
	Reconstituted as per CofA	4°C	1 week
Avoid repeated freeze-thaw cycles.			

www.ptglab.com

Document #: FR-QA118-101 Rev 0
Data Sheet Version #: 1

Proteintech Group, Inc.
5500 Pearl Street, Suite 400 Rosemont, IL 60612
t: 1-888-478-4522; f: 1-312-455-8408
Email: proteintech@ptglab.com