

HumanKine[®] MCP-1 (Recombinant Human)



Animal Component-Free

Human cell expressed

Tag-Free

Endotoxin Free

Product Description

N/A


Alternative Names MCP-1, C-C motif chemokine ligand 2, Monocyte Secretory Protein 1E, CCL2, GDCF-2, HC11, HSMCR30, MCAF, SCYA2, SMC-CF

Source Human Embryonic Kidney cells (HEK293). HEK293-derived MCP-1 protein

Specifications


Test	Method	Specification
Activity	Dose-dependent stimulation of the proliferation of the THP-1 monocyte cell line.	150-750 ng/mL
Molecular Mass	SDS-PAGE	12-14 kDa reduced and non-reduced, monomer, glycosylated
Purity	SDS-PAGE	>95%
Endotoxin	LAL	<1 EU/μg

Activity Data



Recombinant human MCP-1 (HZ-1334) stimulates dose-dependent proliferation of THP-1 human monocyte cell line. Cell number was quantitatively assessed by PrestoBlue[®] Cell Viability Reagent. THP-1 cells were treated with increasing concentrations of recombinant MCP-1 for 144 hours. The EC₅₀ was determined using a 4-parameter non-linear regression model. The EC₅₀

SDS-PAGE



Preparation	
Shipping Temperature	ambient temperature
Formulation	1x PBS pH 7.4, See Certificate of Analysis for details
Reconstitution	Briefly centrifuge the vial before opening. It is recommended to reconstitute the protein to 0.2 mg/mL in sterile 1x PBS pH 7.4 containing 0.1% endotoxin-free recombinant human serum albumin (HSA). Gently swirl or tap vial to mix.

Stability and Storage	Product Form	Temperature Conditions	Storage Time (From Date of Receipt)
	Lyophilized	-20°C to -80°C	Until Expiry Date
	Lyophilized	Room Temperature	2 weeks
	Reconstituted as per CofA	-20°C to -80°C	6 months
	Reconstituted as per CofA	4°C	1 week
Avoid repeated freeze-thaw cycles.			

www.ptglab.com

Document #: FR-QA118-101 Rev 0
Data Sheet Version #:

Proteintech Group, Inc.
5500 Pearl Street, Suite 400 Rosemont, IL 60612
t: 1-888-478-4522; f: 1-312-455-8408
Email: proteintech@ptglab.com