

For Research Use Only

anti-IL6 recombinant VHH, for 1x Cys conjugation



www.ptgcn.com

Catalog Number: IL6Cys1

Basic Information

Catalog Number:
IL6Cys1

Applications:
Conjugation

Host:
Alpaca

Conjugate:
Unconjugated

Type:
Nanobody

Class:
Recombinant

RRID:
AB_3665393

Molecular Weight:
13.7 kDa

Description

IL6Cys1 is an unconjugated recombinant anti IL6 Nanobody (VHH). Suitable for for cysteine conjugation with thiol-reactive reagents, e.g. maleimides. Note: unconjugated VHHs are not suited for usage without prior labeling, since they contain reactive Cysteines. Shipment and storage buffers contain TCEP to keep Cysteines reduced.

Affinity

3 nM

Background

Interleukin-6 (IL-6) is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. IL-6 protein is secreted by a variety of cell types including T cells and macrophages as phosphorylated and variably glycosylated molecule. IL-6 plays an essential role in the final differentiation of B-cells into Ig-secreting cells involved in lymphocyte and monocyte differentiation. It induces myeloma and plasmacytoma growth and induces nerve cells differentiation acts on B-cells, T-cells, hepatocytes, hematopoietic progenitor cells and cells of the CNS. IL-6 is also considered a myokine, a cytokine produced from muscle, and is elevated in response to muscle contraction. IL-6 has been shown to interact with interleukin-6 receptor and glycoprotein 130. Additionally, IL-6 is involved in hematopoiesis, bone metabolism, and cancer progression, and has been defined an essential role in directing transition from innate to acquired immunity.

Storage

Storage:
Store at -20°C

Storage Buffer:
10 mM HEPES pH 7.0, 500 mM NaCl, 0.09% sodium azide


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



BLI analysis of the interaction between varying concentrations of anti-IL-6 VHH and immobilised, biotinylated HumanKine IL-6 (HZ-1019). Fit indicates fitting of data to a 1:1 binding model.