

IHC*easy* ADRP Ready-To-Use IHC Kit

Catalog Number: **KHC0229**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

ADRP (adipocyte differentiation related protein) also known as ADFP, adipophilin, or perilipin-2, is a member of PAT family which is responsible for the transportation of lipids and the formation of lipid droplets. ADRP is localized on the surface of lipid droplets in a variety of tissues and cell lines. ADRP is not detected in undifferentiated cells but increases rapidly to high levels when adipocyte precursors differentiate into adipocytes. Anti-ADRP antibody is a reliable and sensitive marker for lipid droplet. Enhanced expression of ADRP is linked to diseases with abnormal lipid storage, including hepatic steatosis, atherosclerosis and diabetes. Immunohistochemistry of ADRP may facilitate histomorphological diagnosis of these diseases.


Synonyms

ADFP, Adipophilin, ADRP, ADRP/Perilipin 2, Perilipin 2, PLIN2


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using KHC0229 (ADRP/Perilipin 2 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human renal cell carcinoma tissue slide using KHC0229 (ADRP IHC Kit).




Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using KHC0229 (ADRP IHC Kit).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using KHC0229 (ADRP IHC Kit).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0229 (ADRP IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat brown adipose tissue slide using KHC0229 (ADRP IHC Kit).