

IHC*easy* MYO5B Ready-To-Use IHC Kit

Catalog Number: **KHC0356**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

MYO5B is a member of the class V of unconventional, dimeric nonfilamentous myosins, which may be involved in vesicular trafficking. Members of the Rab family of small GTPases, consisting of Rab11a, Rab11b, and Rab25, have recently emerged as potential mediators of vesicle transport by MYO5B. MYO5B is required for polarization of hepatocytes and mutations of MYO5B could disrupt epithelial cell polarity. Loss-of-Function of MYO5B caused by mutation is the main cause of microvillus inclusion disease (MVID). Recently it was reported that the expression of MYO5B was downregulated in gastric cancer and inactivation of MYO5B may contribute to tumorigenesis.


Synonyms

KIAA1119, MYO5B, myosin VB


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using KHC0356 (MYO5B IHC Kit).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0356 (MYO5B IHC Kit).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using KHC0356 (MYO5B IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using KHC0356 (MYO5B IHC Kit).