

IHC*easy* HSPB6 Ready-To-Use IHC Kit

Catalog Number: **KHC0668**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

HSPB6 is a member of the small heat shock superfamily of proteins that function as chaperones to prevent protein misfolding through adenosine triphosphate-independent processes. It is widely expressed but is most abundant in cardiac, skeletal, and smooth muscles. Phosphorylation of HSP20 at Ser16 was elevated in human failing hearts as well as in murine hearts after I/R injury indicating its cardioprotective function.


Synonyms

Heat shock protein beta 6, Hsp20, HSPB6


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC0668 (HSPB6 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using KHC0668 (HSPB6 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using KHC0668 (HSPB6 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse heart tissue slide using KHC0668 (HSPB6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using KHC0668 (HSPB6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat heart tissue slide using KHC0668 (HSPB6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat skeletal muscle tissue slide using KHC0668 (HSPB6 IHC Kit).