

IHCeasy CD117/c-kit Ready-To-Use IHC Kit

Catalog Number: **KHC0818**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

CD117, also named as SCFR, c-Kit and KIT, is a transmembrane tyrosine kinase encoded by the c-Kit proto oncogene. It is a type 3 transmembrane receptor for MGF (mast cell growth factor, also known as stem cell factor). CD117 acts to regulate a variety of biological responses including cell proliferation, apoptosis, chemotaxis and adhesion. Ligand(SCF) binding to the extracellular domain leads to autophosphorylation on several tyrosine residues within the cytoplasmic domain, and activation. Phosphorylation at tyrosine 721 of CD117 allows binding and activation of PI3 kinase. Loss of expression of CD117 appears to be associated with progression of some tumors (melanoma) and autocrine/paracrine stimulation of the kit/SCF system may participate in human solid tumors such as lung, breast, testicular and gynecological malignancies. Mutations in Kit have been found to be important for tumor growth and progression in a variety of cancers including mast cell diseases, gastrointestinal stromal tumor, acute myeloid leukemia, Ewing sarcoma and lung cancer.


Synonyms

KIT/CD117, KIT, CD117, C Kit, EC:2.7.10.1

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung tissue slide using KHC0818 (CD117/c-kit IHC Kit).



Immunohistochemical analysis of paraffin-embedded human appendicitis tissue slide using KHC0818 (CD117/c-kit IHC Kit).