

IHC*easy* PDXK Ready-To-Use IHC Kit

Catalog Number: **KHC0864**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

PDXK encodes for a pyridoxal kinase, which converts inactive B6 vitamers to the active cofactor pyridoxal 5'-phosphate (PLP).


Synonyms

C21orf124, C21orf97, PDXK, PKH, PNK, PRED79, Pyridoxal kinase, Pyridoxine kinase


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC0864 (PDXX IHC Kit).




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC0864 (PDXX IHC Kit).



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using KHC0864 (PDXX IHC Kit).



Immunohistochemical analysis of paraffin-embedded human urothelial carcinoma tissue slide using KHC0864 (PDXX IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC0864 (PDXX IHC Kit).