

IHCeasy[®] TRIM27 Ready-To-Use IHC Kit

Catalog Number: **KHC1532**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Synonyms

Ret finger protein, RFP, RING finger protein 76, RNF76, TRIM27, tripartite motif containing 27, Zinc finger protein RFP


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using KHC1532 (TRIM27 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC1532 (TRIM27 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using KHC1532 (TRIM27 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC1532 (TRIM27 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC1532 (TRIM27 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat colon tissue slide using KHC1532 (TRIM27 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC1532 (TRIM27 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC1532 (TRIM27 IHC Kit).