

# IHC*easy* ZNF24 Ready-To-Use IHC Kit

Catalog Number: **KHC1552**

## General Information

Sample type:  
FFPE tissue  
Cited sample type:  
Reactivity:  
Human, Rat  
Cited Reactivity:

Assay type:  
Immunohistochemistry  
Primary antibody type:  
Rabbit Polyclonal  
Secondary antibody type:  
Polymer-HRP-Goat anti-Rabbit

## Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

## Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


## Background

Zinc finger protein 24 (ZNF24; also known as KOX17, Zfp191 or ZNF191) is a member of the Kruppel-like zinc finger (ZF) transcription factor family, which plays a critical role in cell proliferation and differentiation. ZNF24 protein can be modified with SUMO.

## Synonyms

KOX17, RSG A, Zfp191, Zinc finger protein 191, zinc finger protein 24, Zinc finger protein KOX17, ZNF191, ZNF24, ZSCAN3

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using KHC1552 (ZNF24 IHC Kit).