

IHC*easy* PSMB6 Ready-To-Use IHC Kit

Catalog Number: **KHC1573**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

PSMB6(Proteasome subunit beta type-6) is also named as LMPY, Y and belongs to the peptidase T1B family. The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. It may also catalyzes basal processing of intracellular antigens. It can be up-regulated in anaplastic thyroid cancer cell lines and down-regulated by IFNG/IFN-gamma.


Synonyms

20S proteasome, DELTA, LMPY, Macropain delta chain, Proteasome delta chain, Proteasome subunit beta type 6, Proteasome subunit Y, PSMB6


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using KHC1573 (PSMB6 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using KHC1573 (PSMB6 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC1573 (PSMB6 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using KHC1573 (PSMB6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using KHC1573 (PSMB6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC1573 (PSMB6 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC1573 (PSMB6 IHC Kit).