

IHC*easy* PAICS Ready-To-Use IHC Kit

Catalog Number: **KHC2207**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

PAICS is also named as ADE2(Multifunctional protein), AIRC, PAIS. It is an important bifunctional enzyme in de novo purine biosynthesis in vertebrate with both 5-aminoimidazole ribonucleotide carboxylase (AIRc) and 4-(N-succinylcarboxamide)-5-aminoimidazole ribonucleotide synthetase (SAICARs) activities. It becomes an attractive target for rational anticancer drug design, since rapidly dividing cancer cells rely heavily on the purine de novo pathway for synthesis of adenine and guanine, whereas normal cells favor the salvage pathway.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2207 (PAICS IHC Kit).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC2207 (PAICS IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC2207 (PAICS IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC2207 (PAICS IHC Kit).