

IHC*easy* RBM10 Ready-To-Use IHC Kit

Catalog Number: **KHC2378**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Recombinant
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

RBM10, also named as RNA-binding protein 10, is a 930 amino acid protein, which is in the extranucleolar nucleoplasm constitutes hundreds of nuclear domains, which dynamically change their structures in a reversible manner. RBM10 may be involved in post-transcriptional processing, most probably in mRNA splicing.


Synonyms

RNA binding protein 10, KIAA0122, GPATCH9, GPATC9, G patch domain-containing protein 9


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2378 (RBM10 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using KHC2378 (RBM10 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using KHC2378 (RBM10 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC2378 (RBM10 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using KHC2378 (RBM10 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC2378 (RBM10 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC2378 (RBM10 IHC Kit).