

# IHC*easy* DUSP2 Ready-To-Use IHC Kit

Catalog Number: **KHC2394**

## General Information

Sample type:  
FFPE tissue  
Cited sample type:  
Reactivity:  
Human, Rat  
Cited Reactivity:

Assay type:  
Immunohistochemistry  
Primary antibody type:  
Rabbit Polyclonal  
Secondary antibody type:  
Polymer-HRP-Goat anti-Rabbit

## Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

## Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


## Background

DUSP2 belongs to the sub-family of nuclear specific DUSPs that includes DUSP1, 4 and 5. DUSP2 is highly expressed in immune cells, in particular, of patients with inflammatory arthritis, such as in rheumatoid arthritis (RA), and Dusp2 knockout mice present a significant reduction in inflammatory responses. It has been demonstrated that DUSP2 is a downstream target of the tumor suppressor p53 that is required for the cellular response to oxidative stress. DUSP2 is also involved in transcription factor E2F1-mediated apoptotic signaling. Like other members of the DUSP family, DUSP2 interacts with MAPKs; However, the target protein for its phosphatase activity is controversial.


## Synonyms

PAC1, PAC 1, EC:3.1.3.16, Dual specificity protein phosphatase PAC-1, dual specificity phosphatase 2

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human urothelial carcinoma tissue slide using KHC2394 (DUSP2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat testis tissue slide using KHC2394 (DUSP2 IHC Kit).