

IHC*easy* PCBP2 Ready-To-Use IHC Kit

Catalog Number: **KHC2486**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

Poly(rC)-binding protein 2 (PCBP2) is a single-stranded nucleic acid binding protein that binds preferentially to oligo dC. It is also named Heterogeneous nuclear ribonucleoprotein E2, binds also poly(rU). It negatively regulates cellular antiviral responses mediated by MAVS signaling. It acts as an adapter between MAVS and the E3 ubiquitin ligase ITCH, therefore triggering MAVS ubiquitination and degradation.


Synonyms

hnRNP E2,PCBP2,Alpha CP2,Alpha-CP2,HNRPE2


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2486 (PCBP2 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human rectal cancer tissue slide using KHC2486 (PCBP2 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC2486 (PCBP2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC2486 (PCBP2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat cerebellum tissue slide using KHC2486 (PCBP2 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC2486 (PCBP2 IHC Kit).