

IHC*easy* FTSJ3 Ready-To-Use IHC Kit

Catalog Number: **KHC2657**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

FTSJ3 is an RNA 2'-O-methyltransferase recruited by HIV to avoid innate immune sensing.


Synonyms

FTSJ3, EC:2.1.1.-, EPCS3, FTSJ 3, pre-rRNA 2'-O-ribose

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2657 (FTSJ3 IHC Kit).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC2657 (FTSJ3 IHC Kit).