

IHCeasy[®] TMEM111 Ready-To-Use IHC Kit

Catalog Number: **KHC2984**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

TMEM111, also names as EMC3, is part of the endoplasmic reticulum-associated secretory pathway. TMEM111 was required for murine pulmonary surfactant synthesis and lung function at birth.


Synonyms

EMC3, POB


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC2984 (TMEM111 IHC Kit).




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2984 (TMEM111 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC2984 (TMEM111 IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using KHC2984 (TMEM111 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC2984 (TMEM111 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat pancreas tissue slide using KHC2984 (TMEM111 IHC Kit).