

IHCeasy[®] UVRAG Ready-To-Use IHC Kit

Catalog Number: KHC2997

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

UVRAG is a versatile protein that is involved in the regulation of different cellular pathways implicated in membrane trafficking, and is also involved in the regulation of the COPI-dependent retrograde transport from Golgi and the endoplasmic reticulum by associating with the NRZ complex; Its function is dependent on its binding to phosphatidylinositol 3-phosphate (PtdIns3P).


Synonyms

DHTX, p63, UV radiation resistance-associated gene protein


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2997 (UVRAG IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC2997 (UVRAG IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC2997 (UVRAG IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC2997 (UVRAG IHC Kit).