

IHCeasy[®] FAM13A Ready-To-Use IHC Kit

Catalog Number: **KHC3099**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

FAM13A, also named as FAM13A1 and KIAA0914, belongs to the FAM13 family. It is associated with chronic obstructive pulmonary disease.


Synonyms

FAM13A1, KIAA0914, Protein FAM13A


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human heart tissue slide using KHC3099 (FAM13A IHC Kit).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using KHC3099 (FAM13A IHC Kit).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC3099 (FAM13A IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using KHC3099 (FAM13A IHC Kit).