

For Research Use Only

ISG15-Trap Magnetic Agarose



www.ptgcn.com

Catalog Number: **MG-3001MA**

Basic Information

Catalog Number:
MG-3001MA

Applications:
IP, Co-IP

Host:
Alpaca

Type:
VHH

Class:
Recombinant - Animal free production

Description

The ChromoTek ISG15-Trap Magnetic Agarose consists of an anti-ISG15 VHH (Nanobody), which is coupled to magnetic agarose beads. It can be used for the immunoprecipitation of ISG15 tagged proteins from cell extracts of various organisms.

Specificity/Target

/

Elution buffer

2x SDS-sample buffer (Lämmli), 200 mM glycine pH 2.5

Affinity

70 nM

Storage

Storage:
Upon arrival store at +4°C / do not freeze!

Storage Buffer:
20% Ethanol


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



ISG15-Trap magnetic agarose (MG-3001MA) was used to immunoprecipitate endogenous ISG15 and ISGylated proteins from untreated human HEPG2 cells (baseline) and interferon- β treated HepG2 cells (induced). Lysis was achieved with RIPA buffer. For each IP, samples of the input lysate (I), non-bound flow-through (FT), and bound (B) fractions were analyzed via Coomassie stained SDS-PAGE & Western blot. Anti-ISG15 polyclonal antibody (15981-1-AP) Multi-rAb® HRP-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR001) were used in the Western blot analysis. The Trap shows efficient IP of endogenous ISG15 and ISGylated proteins with low background.

ISG15-Trap magnetic agarose (MG-3001MA) was used to immunoprecipitate endogenous ISG15 and ISGylated proteins from untreated human HEPG2 cells (baseline) and interferon- β treated HepG2 cells (induced). Lysis was achieved with standard Lysis buffer. For each IP, samples of the input lysate (I), non-bound flow-through (FT), and bound (B) fractions were analyzed via Coomassie stained SDS-PAGE & Western blot. Anti-ISG15 polyclonal antibody (15981-1-AP) and Multi-rAb® HRP-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR001) were used in the Western blot analysis. The Trap shows efficient IP of endogenous ISG15 and ISGylated proteins with low background.