

ICAM-1/CD54 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00028-2**

Capture Antibody Information

Catalog Number:
83069-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230342F1

Reactivity:
human

GenBank:
BC015969

Conjugate:
Unconjugated

Full name:
intercellular adhesion molecule 1

Gene ID:
3383

Detection Antibody Information

Catalog Number:
83069-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230342H4

Reactivity:
human

GenBank:
BC015969

Conjugate:
Unconjugated

Full name:
intercellular adhesion molecule 1

Gene ID:
3383

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00028-2 targets ICAM-1/CD54 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ICAM-1/CD54 Recombinant antibody, PBS Only (Capture) 83069-1-PBS (230342F1). 100 µg. Concentration 1 mg/mL.

Detection antibody: ICAM-1/CD54 Recombinant antibody, PBS Only (Detector) 83069-4-PBS (230342H4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00028-2, ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83069-1-PBS. Detection antibody: 83069-4-PBS. Standard: Eg0570. Range: 0.625-20 ng/mL