

# IFT140 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00029-2

## Capture Antibody Information

<b>Catalog Number:</b> 82947-4-PBS	<b>Clone ID:</b> 230231A6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> intraflagellar transport 140 homolog (Chlamydomonas)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC035577	<b>Gene ID:</b> 9742
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag11119	

## Detection Antibody Information

<b>Catalog Number:</b> 82947-6-PBS	<b>Clone ID:</b> 230231H2	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> intraflagellar transport 140 homolog (Chlamydomonas)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC035577	<b>Gene ID:</b> 9742
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag11119	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 2.5-20 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00029-2 targets IFT140 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IFT140 Recombinant antibody, PBS Only (Capture) 82947-4-PBS (230231A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: IFT140 Recombinant antibody, PBS Only (Detector) 82947-6-PBS (230231H2). 100 µg. Concentration 1 mg/mL.

Alternative IFT140 matched antibody pairs: MP00029-1, MP00029-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00029-2, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-4-PBS. Detection antibody: 82947-6-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL.