

For Research Use Only

IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00035-3**

Capture Antibody Information

Catalog Number:
83008-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240008G6
Reactivity:
human
GenBank:
BC067498

Conjugate:
Unconjugated
Full name:
interleukin 12B (natural killer cell
stimulatory factor 2, cytotoxic
lymphocyte maturation factor 2, p40)
Gene ID:
3593

Detection Antibody Information

Catalog Number:
83008-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240008F4
Reactivity:
human
GenBank:
BC067498

Conjugate:
Unconjugated
Full name:
interleukin 12B (natural killer cell
stimulatory factor 2, cytotoxic
lymphocyte maturation factor 2, p40)
Gene ID:
3593

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00035-3 targets IL-12/IL-23 p40 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-12/IL-23 p40 Recombinant antibody, PBS Only (Capture) 83008-3-PBS (240008G6). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-12/IL-23 p40 Recombinant antibody, PBS Only (Detector) 83008-1-PBS (240008F4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00035-3, IL12B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83008-3-PBS. Detection antibody: 83008-1-PBS. Standard: Eg0431. Range: 0.625-20 ng/mL