

CAMP Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00042-1

Capture Antibody Information

Catalog Number:
82940-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230157C4
Reactivity:
Human
GenBank:
BC055089
Immunogen Catalog Number:
Ag2622

Conjugate:
Unconjugated
Full name:
cathelicidin antimicrobial peptide
Gene ID:
820

Detection Antibody Information

Catalog Number:
82940-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230157E11
Reactivity:
Human
GenBank:
BC055089
Immunogen Catalog Number:
Ag2622

Conjugate:
Unconjugated
Full name:
cathelicidin antimicrobial peptide
Gene ID:
820

Applications

Tested Applications:
Sandwich ELISA

Range:
0.781-50ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00042-1 targets CAMP in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: CAMP Recombinant antibody, PBS Only (Capture) 82940-3-PBS (230157C4). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CAMP Recombinant antibody, PBS Only (Detector) 82940-5-PBS (230157E11). 100 μ g. Concentration 1 mg/mL.

Alternative CAMP matched antibody pairs: MP00042-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve ofMP00042-1, CAMP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82940-3-PBS. Detection antibody: 82940-5-PBS. Standard: Ag2622. Range: 0.781-50ng/mL.