

# CTAGE5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00052-1

## Capture Antibody Information

**Catalog Number:**  
83017-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230262E6  
**Reactivity:**  
Human  
**GenBank:**  
NM\_005930

**Conjugate:**  
Unconjugated  
**Full name:**  
CTAGE family, member 5  
**Gene ID:**  
4253

## Detection Antibody Information

**Catalog Number:**  
83017-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230262E12  
**Reactivity:**  
Human  
**GenBank:**  
NM\_005930

**Conjugate:**  
Unconjugated  
**Full name:**  
CTAGE family, member 5  
**Gene ID:**  
4253

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.156-10 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00052-1 targets CTAGE5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CTAGE5 Recombinant antibody, PBS Only (Capture) 83017-1-PBS (230262E6). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: CTAGE5 Recombinant antibody, PBS Only (Detector) 83017-3-PBS (230262E12). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative CTAGE5 matched antibody pairs: MP00052-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00052-1, CTAGE5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83017-1-PBS. Detection antibody: 83017-3-PBS. Standard: SY01342. Range: 0.156-10 ng/mL.