

TMEM132E Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00054-1

Capture Antibody Information

Catalog Number:
82983-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230525D9
Reactivity:
Human
GenBank:
NM_001304438
Immunogen Catalog Number:
Ag26772

Conjugate:
Unconjugated
Full name:
transmembrane protein 132E
Gene ID:
124842

Detection Antibody Information

Catalog Number:
82983-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230525F7
Reactivity:
Human
GenBank:
NM_001304438
Immunogen Catalog Number:
Ag26772

Conjugate:
Unconjugated
Full name:
transmembrane protein 132E
Gene ID:
124842

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00054-1 targets TMEM132E in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMEM132E Recombinant antibody, PBS Only (Capture) 82983-5-PBS (230525D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMEM132E Recombinant antibody, PBS Only (Detector) 82983-3-PBS (230525F7). 100 µg. Concentration 1 mg/mL.

Alternative TMEM132E matched antibody pairs: MP00054-2, MP00054-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00054-1, TMEM132E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82983-5-PBS. Detection antibody: 82983-3-PBS. Standard: Ag26772. Range: 0.625-20 ng/mL.