

FANCM Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00057-1

Capture Antibody Information

Catalog Number:
82897-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230233G9
Reactivity:
Human
GenBank:
BC140776
Immunogen Catalog Number:
Ag34447

Conjugate:
Unconjugated
Full name:
Fanconi anemia, complementation group M
Gene ID:
57697

Detection Antibody Information

Catalog Number:
82897-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230233D10
Reactivity:
Human
GenBank:
BC140776
Immunogen Catalog Number:
Ag34447

Conjugate:
Unconjugated
Full name:
Fanconi anemia, complementation group M
Gene ID:
57697

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00057-1 targets FANCM in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FANCM Recombinant antibody, PBS Only (Capture) 82897-6-PBS (230233G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: FANCM Recombinant antibody, PBS Only (Detector) 82897-5-PBS (230233D10). 100 µg. Concentration 1 mg/mL.

Alternative FANCM matched antibody pairs: MP00057-2, MP00057-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00057-1, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-6-PBS. Detection antibody: 82897-5-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL