

For Research Use Only

# EPOR Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00065-1

## Capture Antibody Information

Catalog Number:  
83056-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230117G3

Reactivity:  
human

GenBank:  
NM\_000121.4

Immunogen Catalog Number:  
Ag33956

Conjugate:  
Unconjugated

Full name:  
erythropoietin receptor

Gene ID:  
2057

## Detection Antibody Information

Catalog Number:  
83056-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
230117F9

Reactivity:  
human

GenBank:  
NM\_000121.4

Immunogen Catalog Number:  
Ag33956

Conjugate:  
Unconjugated

Full name:  
erythropoietin receptor

Gene ID:  
2057

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
0.625-40 ng/mL (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00065-1 targets EPOR in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: EPOR Recombinant antibody, PBS Only (Capture) 83056-4-PBS (230117G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: EPOR Recombinant antibody, PBS Only (Detector) 83056-3-PBS (230117F9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00065-1, EPOR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83056-4-PBS. Detection antibody: 83056-3-PBS. Standard: Ag33956. Range: 0.625-40 ng/mL