

STT3A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00080-1**

Capture Antibody Information

Catalog Number:
83077-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230426G7
Reactivity:
Human
GenBank:
BC020965
Immunogen Catalog Number:
Ag31994

Conjugate:
Unconjugated
Full name:
STT3, subunit of the oligosaccharyltransferase complex, homolog A (S. cerevisiae)
Gene ID:
3703

Detection Antibody Information

Catalog Number:
83077-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230426D2
Reactivity:
Human
GenBank:
BC020965
Immunogen Catalog Number:
Ag31994

Conjugate:
Unconjugated
Full name:
STT3, subunit of the oligosaccharyltransferase complex, homolog A (S. cerevisiae)
Gene ID:
3703

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00080-1 targets STT3A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STT3A Recombinant antibody, PBS Only (Capture) 83077-3-PBS (230426G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: STT3A Recombinant antibody, PBS Only (Detector) 83077-2-PBS (230426D2). 100 µg. Concentration 1 mg/mL.

Alternative STT3A matched antibody pairs: MP00080-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00080-1, STT3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83077-3-PBS. Detection antibody: 83077-2-PBS. Standard: Ag31994. Range: 0.625-20 ng/mL