

ACVR1B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00098-2

Capture Antibody Information

Catalog Number: 83026-1-PBS	Clone ID: 230260C5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: activin A receptor, type IB
Isotype: IgG	GenBank: BC000254	Gene ID: 91
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0129	

Detection Antibody Information

Catalog Number: 83026-4-PBS	Clone ID: 230260H5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: activin A receptor, type IB
Isotype: IgG	GenBank: BC000254	Gene ID: 91
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0129	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00098-2 targets ACVR1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ACVR1B Recombinant antibody, PBS Only (Capture) 83026-1-PBS (230260C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: ACVR1B Recombinant antibody, PBS Only (Capture/Detector) 83026-4-PBS (230260H5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00098-2, ACVR1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83026-1-PBS. Detection antibody: 83026-4-PBS. Standard: Ag0129. Range: 0.313-20 ng/mL