

ARID1B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00103-3**

Capture Antibody Information

Catalog Number:
82979-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230222H3
Reactivity:
Human
GenBank:
NM_017519
Immunogen Catalog Number:
Ag30643

Conjugate:
Unconjugated
Full name:
AT rich interactive domain 1B (SWI1-like)
Gene ID:
57492

Detection Antibody Information

Catalog Number:
82979-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230222B7
Reactivity:
Human
GenBank:
NM_017519
Immunogen Catalog Number:
Ag30643

Conjugate:
Unconjugated
Full name:
AT rich interactive domain 1B (SWI1-like)
Gene ID:
57492

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00103-3 targets ARID1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ARID1B Recombinant antibody, PBS Only (Capture) 82979-2-PBS (230222H3). 100 µg. Concentration 1 mg/mL.

Detection antibody: ARID1B Recombinant antibody, PBS Only (Detector) 82979-3-PBS (230222B7). 100 µg. Concentration 1 mg/mL.

Alternative ARID1B matched antibody pairs: MP00103-1, MP00103-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00103-3, ARID1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82979-2-PBS. Detection antibody: 82979-3-PBS. Standard: Ag30643. Range: 0.625-20 ng/mL