

# STK39 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00106-3**

## Capture Antibody Information

<b>Catalog Number:</b> 83081-2-PBS	<b>Clone ID:</b> 230382G9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human, mouse, rat	<b>Full name:</b> serine threonine kinase 39 (STE20/SPS1 homolog, yeast)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC166614	<b>Gene ID:</b> 27347
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag34503	

## Detection Antibody Information

<b>Catalog Number:</b> 83081-1-PBS	<b>Clone ID:</b> 230382C5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> serine threonine kinase 39 (STE20/SPS1 homolog, yeast)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC166614	<b>Gene ID:</b> 27347
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag34503	

## Applications

<b>Tested Applications:</b> Sandwich ELISA	<b>Range:</b> 31.3-2000 pg/mL (Sandwich ELISA)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	---

## Product Information

MP00106-3 targets STK39 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: STK39 Recombinant antibody, PBS Only (Capture) 83081-2-PBS (230382G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: STK39 Recombinant antibody, PBS Only (Detector) 83081-1-PBS (230382C5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Sandwich ELISA standard curve of MP00106-3, Human STK39 Recombinant Matched Antibody Pair - PBS only. 83081-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34503. 83081-1-PBS was HRP conjugated as the detection antibody. Range: 31.3-2000 pg/mL.



The mean STK39 concentration was determined to be 2,623.5 pg/mL in HepG2 cell extract based on a 4.2 mg/mL extract load and 9,245.2 pg/mL in Jurkat cell extract based on a 4.6 mg/mL extract load.