

Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00116-1**

Capture Antibody Information

Catalog Number:
83064-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230336H6
Reactivity:
mouse
GenBank:
NM_134248
Immunogen Catalog Number:
Eg0571

Conjugate:
Unconjugated
Full name:
hepatitis A virus cellular receptor 1
Gene ID:
171283

Detection Antibody Information

Catalog Number:
83064-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230336D2
Reactivity:
mouse
GenBank:
NM_134248
Immunogen Catalog Number:
Eg0571

Conjugate:
Unconjugated
Full name:
hepatitis A virus cellular receptor 1
Gene ID:
171283

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00116-1 targets KIM-1/HAVCR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only 83064-6-PBS (230336H6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only 83064-3-PBS (230336D2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00116-1, Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83064-6-PBS. Detection antibody: 83064-3-PBS. Standard: Eg0571. Range: 0.625-10 ng/mL