

TMEM20 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00120-3

Capture Antibody Information

Catalog Number: 83045-3-PBS	Clone ID: 230298E1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: transmembrane protein 20
Isotype: IgG	GenBank: BC104814	Gene ID: 159371
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34357	

Detection Antibody Information

Catalog Number: 83045-1-PBS	Clone ID: 230298C11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: transmembrane protein 20
Isotype: IgG	GenBank: BC104814	Gene ID: 159371
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34357	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00120-3 targets TMEM20 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMEM20 Recombinant antibody, PBS Only (Capture) 83045-3-PBS (230298E1). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMEM20 Recombinant antibody, PBS Only (Detector) 83045-1-PBS (230298C11). 100 µg. Concentration 1 mg/mL.

Alternative TMEM20 matched antibody pairs: MP00120-1, MP00120-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00120-3, TMEM20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83045-3-PBS. Detection antibody: 83045-1-PBS. Standard: Ag34357. Range: 0.625-20 ng/mL.