

SLC35F4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00122-1

Capture Antibody Information

Catalog Number:
83085-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230360G9
Reactivity:
Human
GenBank:
NM_001206920
Immunogen Catalog Number:
Ag34518

Conjugate:
Unconjugated
Full name:
solute carrier family 35, member F4
Gene ID:
341880

Detection Antibody Information

Catalog Number:
83085-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230360H1
Reactivity:
Human
GenBank:
NM_001206920
Immunogen Catalog Number:
Ag34518

Conjugate:
Unconjugated
Full name:
solute carrier family 35, member F4
Gene ID:
341880

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00122-1 targets SLC35F4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC35F4 Recombinant antibody, PBS Only (Capture) 83085-1-PBS (230360G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC35F4 Recombinant antibody, PBS Only (Detector) 83085-2-PBS (230360H1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00122-1, SLC35F4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83085-1-PBS. Detection antibody: 83085-2-PBS. Standard: Ag34518. Range: 1.25-20 ng/mL