

NFX1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00127-3**

Capture Antibody Information

Catalog Number:
82924-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230211G5
Reactivity:
human
GenBank:
BC012151
Immunogen Catalog Number:
Ag9448

Conjugate:
Unconjugated
Full name:
nuclear transcription factor, X-box binding 1
Gene ID:
4799

Detection Antibody Information

Catalog Number:
82924-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230211H6
Reactivity:
human
GenBank:
BC012151
Immunogen Catalog Number:
Ag9448

Conjugate:
Unconjugated
Full name:
nuclear transcription factor, X-box binding 1
Gene ID:
4799

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00127-3 targets NFX1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFX1 Recombinant antibody, PBS Only 82924-5-PBS (230211G5). 100 µg. Concentration 1 mg/mL.

Detection antibody: NFX1 Recombinant antibody, PBS Only 82924-6-PBS (230211H6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00127-3, NFX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82924-5-PBS. Detection antibody: 82924-6-PBS. Standard: Ag9448. Range: 0.625-20 ng/mL