

uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00129-1

Capture Antibody Information

Catalog Number: 82953-2-PBS	Clone ID: 230148B10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, mouse	Full name: PLAUR
Isotype: IgG	GenBank: BC002788	Gene ID: 5329
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28056	

Detection Antibody Information

Catalog Number: 82953-3-PBS	Clone ID: 230148E10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: PLAUR
Isotype: IgG	GenBank: BC002788	Gene ID: 5329
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28056	

Applications

Tested Applications: Sandwich ELISA	Range: 125-8000 pg/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	--

Product Information

MP00129-1 targets uPAR/CD87 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: uPAR/CD87 Recombinant antibody, PBS Only (Capture/Detector) 82953-2-PBS (230148B10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: uPAR/CD87 Recombinant antibody, PBS Only (Detector) 82953-3-PBS (230148E10). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP00129-1, uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82953-2-PBS. Detection antibody: 82953-3-PBS. Standard: Ag28056. Range: 125-8000 pg/mL