

Cytokeratin 1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00134-1

Capture Antibody Information

Catalog Number:

82972-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230261E1

Reactivity:

human

GenBank:

NM_006121

Conjugate:

Unconjugated

Full name:

keratin 1

Gene ID:

3848

Detection Antibody Information

Catalog Number:

82972-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230261D12

Reactivity:

human

GenBank:

NM_006121

Conjugate:

Unconjugated

Full name:

keratin 1

Gene ID:

3848

Applications

Tested Applications:

Sandwich ELISA

Range:

0.313-20 ng/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00134-1 targets Cytokeratin 1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Cytokeratin 1-specific Recombinant antibody, PBS Only (Capture) 82972-3-PBS (230261E1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Cytokeratin 1-specific Recombinant antibody, PBS Only (Detector) 82972-2-PBS (230261D12). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP00134-1, Cytokeratin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82972-3-PBS. Detection antibody: 82972-2-PBS. Standard: SY00109. Range: 0.313-20 ng/mL.