

SLC25A28 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00142-1**

Capture Antibody Information

Catalog Number:
83127-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230392D4

Reactivity:
Human

GenBank:
BC076399

Immunogen Catalog Number:
Ag34507

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 28

Gene ID:
81894

Detection Antibody Information

Catalog Number:
83127-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230392C1

Reactivity:
Human

GenBank:
BC076399

Immunogen Catalog Number:
Ag34507

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 28

Gene ID:
81894

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00142-1 targets SLC25A28 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A28 Recombinant antibody, PBS Only (Capture) 83127-4-PBS (230392D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC25A28 Recombinant antibody, PBS Only (Capture/Detector) 83127-2-PBS (230392C1). 100 µg. Concentration 1 mg/mL.

Alternative SLC25A28 matched antibody pairs: MP00142-2, MP00142-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00142-1, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-4-PBS. Detection antibody: 83127-2-PBS. Standard: Ag34507. Range: 0.625-20 ng/mL.