

PGF/PIGF Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00146-1

Capture Antibody Information

Catalog Number:
83170-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230532H3
Reactivity:
human
GenBank:
BC007789
Immunogen Catalog Number:
Ag0969

Conjugate:
Unconjugated
Full name:
placental growth factor
Gene ID:
5228

Detection Antibody Information

Catalog Number:
83170-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230532C12
Reactivity:
human
GenBank:
BC007789
Immunogen Catalog Number:
Ag0969

Conjugate:
Unconjugated
Full name:
placental growth factor
Gene ID:
5228

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00146-1 targets PGF/PIGF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PGF/PIGF Recombinant antibody, PBS Only (Capture) 83170-4-PBS (230532H3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PGF/PIGF Recombinant antibody, PBS Only (Capture/Detector) 83170-2-PBS (230532C12). 100 μ g. Concentration 1 mg/mL.

Alternative PGF/PIGF matched antibody pairs: MP00146-2, MP00146-3, MP50063-1, MP50063-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00146-1, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-4-PBS. Detection antibody: 83170-2-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL