

PGF/PIGF Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00146-3

Capture Antibody Information

Catalog Number: 83170-2-PBS	Clone ID: 230532C12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: placental growth factor
Isotype: IgG	GenBank: BC007789	Gene ID: 5228
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0969	

Detection Antibody Information

Catalog Number: 83170-1-PBS	Clone ID: 230532A8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: placental growth factor
Isotype: IgG	GenBank: BC007789	Gene ID: 5228
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0969	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
-----------------------------------------------	--------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP00146-3 targets PGF/PIGF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PGF/PIGF Recombinant antibody, PBS Only (Capture/Detector) 83170-2-PBS (230532C12). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PGF/PIGF Recombinant antibody, PBS Only (Detector) 83170-1-PBS (230532A8). 100 μ g. Concentration 1 mg/mL.

Alternative PGF/PIGF matched antibody pairs: MP00146-1, MP00146-2, MP50063-1, MP50063-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00146-3, PGF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83170-2-PBS. Detection antibody: 83170-1-PBS. Standard: Ag0969. Range: 0.625-20 ng/mL