

Nectin-4/PVRL4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00147-1

Capture Antibody Information

Catalog Number:
83223-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240009E4
Reactivity:
human
GenBank:
BC010423

Conjugate:
Unconjugated
Full name:
poliovirus receptor-related 4
Gene ID:
81607
Affinity:
 $K_D = 8.49 \times 10^{-10}$
 $K_{off} = 2.20 \times 10^{-4}$
 $K_{on} = 2.59 \times 10^5$

Detection Antibody Information

Catalog Number:
83223-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240009F2
Reactivity:
human
GenBank:
BC010423

Conjugate:
Unconjugated
Full name:
poliovirus receptor-related 4
Gene ID:
81607
Affinity:
 $K_D = 2.28 \times 10^{-10}$
 $K_{off} = 2.98 \times 10^{-5}$
 $K_{on} = 2.62 \times 10^5$

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00147-1 targets Nectin-4/PVRL4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: nectin 4 Recombinant antibody, PBS Only (Capture) 83223-1-PBS (240009E4). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Nectin-4/PVRL4 Recombinant antibody, PBS Only (Detector) 83223-4-PBS (240009F2). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00147-1, Nectin-4/PVRL4 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83223-1-PBS. Detection antibody: 83223-4-PBS. Standard: Eg0415. Range: 0.625-10 ng/mL.